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Abstract

Continuous-flow (on-line) system has revolutionized the study of natural stable isotope variations in organic matter and organic compounds. In particular, carbon and nitrogen isotope analysis of bulk samples using elemental analyzer/isotope ratio mass spectrometer (EA/IRMS) has been widely used in various field of studies, which allows a simple and rapid analysis of the isotope ratios for all types of organic and a wide range of inorganic samples, with a small amount of sample materials compared to the traditional off-line Dual-Inlet method. However, stability, accuracy, and precision on the observed isotope ratios strongly reflect various factors associated with continuous-flow system, such as peak intensity (*i.e.*, sample weight), memory effect, and He dilution ratio. Unfortunately, the complete removal of such uncertainties should not be realistic, particularly in routine isotope measurements of a great number of samples. Here, we review the principles of stable isotope analysis on the traditional off-line Dual-Inlet and continuous-flow EA/IRMS methods, and also show an applicable calibration and standardization sequence. We hope that this paper is useful for routine application of EA/IRMS technique to various studies in geochemistry, ecology, biology, food science, and other area of sciences.

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*Carbon and Nitrogen stable isotope analysis by EA/IRMS

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Yaeko Suzuki: Japan Certification Services, Inc. 4-5-17 Chigasaki-Higashi, Tsuzuki, Yokohama, Kanagawa, 224-0033, Japan

Corresponding author. Rie Sato

e-mail: r-sato@shoko.co.jp Tel: +81-480-35-2915; Fax: +81-480-35-2930

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Rie Sato: SI Science CO., LTD. 473-3 Hongo, Sugito-Machi, Kitakatsushika-Gun, Saitama, 345-0023, Japan